



Biopex A Bioreactor | Fermenter
Technical Data Sheet

Biopex A Bioreactor | Fermenter is a laboratory equipment for microbial fermentations and cell culture applications.

Unopex offers different configurations:

- **Biopex A**
Single configuration with one culture vessel and one control tower.
- **Biopex A twin**
Twin configuration with two culture vessels and one control tower.
- **Biopex A poly**
Poly configuration with up to 6 culture vessels, one supply tower per vessel and one control tower.



Standard Instrument

Code	Instrument
A50122	Biopex A
A52145	Biopex A twin
A55450	Biopex A poly

Instruments are delivered ready to use.

Configuration

Component	Biopex A	Biopex A twin	Biopex A poly
Operation with one culture vessel	✓		
Independent operation of two culture vessels		✓	
Independent operation of up to six culture vessels			✓
Culture vessel options single-wall or double-wall 1 L 2 L 5 L 10 L	✓	✓	✓
Control tower with touch screen operating panel	✓	✓	✓
User-friendly software for monitoring and operating the entire process	✓	✓	✓
Supply tower			✓
pH sensor	✓	✓	✓
pO ₂ sensor	✓	✓	✓
Level/foam sensor	✓	✓	✓
Peristaltic pump	✓	✓	✓
Stirrer with servo motor	✓	✓	✓
6-blade Rushton impeller	✓	✓	✓
3-blade segment impeller	✓	✓	✓
Exhaust cooler	✓	✓	✓
Inlet/exhaust gas filters	✓	✓	✓
Triple/quadruple port	✓	✓	✓
Heating blanket	✓	✓	✓
Cooling finger	✓	✓	✓
Mass flow controllers	✓	✓	✓
Ring sparger	✓	✓	✓
Micro sparger	✓	✓	✓
Removable baffle	✓	✓	✓
Reagent bottles	✓	✓	✓
Vessel stand	✓	✓	✓
Adjustable dip tube	✓	✓	✓
Resterilizable sampler	✓	✓	✓
Operation & Maintenance Manual	✓	✓	✓

Scope of supply might change according to specific offers, microbial or cell culture applications.

Principle of Operation

Fully-automated touch screen operating panel with smart process control software ensures Biopex A bioreactor | fermenter to be operated in the required process conditions.

- Stirrer operates with the energy it receives from the motor.
- Gassing is carried out by sparger into the culture medium;
 - individual gas flow paths
 - mix of gas types
- pH is controlled by acid and base pumps.
- pO₂ is controlled by multiple cascade system (air, O₂, N₂, stirrer speed, feed).
- Level/foam is controlled by peristaltic pump AFoam/Level.
- Feed is controlled by Feed/Harvest pump.
- Temperature is controlled by heating blanket & cooling finger (for single-wall vessel) or thermal circulator (for double-wall vessel)

Technical Data

Working volume	1 L	0,3 1 L
	2 L	0,6 2 L
	5 L	0,8 5 L
	10 L	1,4 10 L
Pumps	Acid, Base, AFoam/Level, Feed/Harvest	
Operating panel	Touchscreen	
Range of speed	1 L	0 2000 rpm
	2 L	0 2000 rpm
	5 L	0 1500 rpm
	10 L	0 1500 rpm
Impeller type	6-blade Rushton impeller for microbial version 3-blade segment impeller for cell culture version	
Gas control	Mass flow controllers with valves for Air, O ₂ for microbial version Air, O ₂ , N ₂ , CO ₂ , head space gassing for cell culture version	
Connections	Quick connect couplings	
Sensors	pH temperature - pO ₂ level-foam	
Temperature control range	20 60 °C for double wall culture vessels (20 80 °C is available as an option)	
Computer connection Data transfer	Ethernet / usb	
Connection voltage	220 V - 50 Hz	
Material of construction for product contact parts	Heat resistant borosilicate glass, Stainless steel AISI 316L, Silicone, PTFE/EPDM	
Culture vessel size in autoclave (W x D x H)	1 L	170 x 190 x 265 mm
	2 L	270 x 210 x 305 mm
	5 L	295 x 245 x 405 mm
	10 L	325 x 310 x 560 mm

Options & Accessories

Component
Integrated thermal circulator for temperature control
Additional mass flow controllers/rotameters
External pumps
Desktop computer & printer
Automatic weight control
Gravimetric flow control
Gravimetric level control
Tangential flow filtration (TFF)
Turbidity sensor
Conductivity sensor
ORP measurement
IQ/OQ documentation